




1. Remove items from box and place on a clean surface. Peel off protective cover from visor.



2. Ensure headset is positioned correctly with "top side" facing up. Align holes in visor with pins in headset (left to right).



3. Attach visor to headset by snap-fitting each prong to hole (left to right).



4. Thread headband band through tabs in headset.




5. Apply protective face shield to head and secure in place with headband.

# Reusable Protective Face Shield

Personal protective equipment (PPE) designed to help protect your face, mouth and eyes from spray and splatter.

## Product Features:

- Fast assembly and easy cleaning
- Adjustable Elastic band offers secure fit



## DO NOT AUTOCLAVE

Each set includes one each of these component parts:



## SAFETY TIPS

- Each headset should be cleaned between uses
- Parts should not be shared between persons without submitting for sterilization

## Headset

### Materials

**Headband** - Black ABS

**Foam Liner** - Closed Cell, Crosslined Polyethylene foam

**Foam Insert** - Closed cell, crosslinked Polyethylene foam

### Usage

Multi-Use

### Sterilization

PACKAGED NON-STERILE

## Sterilization/Cleaning Instructions

Bleach (sodium hypochlorite)  
0.5% (1:10 bleach: water), 2 mins  
Soap/detergent + water  
Cleaning of dirt only (not sterilizing)  
air dry

## Please Recycle



PET  
VISOR



PP  
HEADSET

## Visor

### Dimensions

240mm x 317 mm  
240mm x 360 mm

### Materials

.51mm Thick PETG

### Usage

Multi-Use

### Sterilization

PACKAGED NON-STERILE

## Sterilization/Cleaning Instructions

Isopropyl alcohol - 75-99% 5 minutes  
Bleach (sodium hypochlorite) 0.5%  
(1:10 bleach: water), 2 mins  
Ethanol 70-80% maximum, 5 mins

## Headband

### Dimensions

Adjustable Elastic Band

### Materials

Polyester Woven Elastic

### Usage

Single-Use

### Sterilization

PACKAGED NON-STERILE

## Cleaning Instructions

Launder Wash / Standard Dry Cycle

## Flexfit

The headband is designed for a adjustable fit. Re-position the headband for a more comfortable fit.

100% Satisfaction Guaranteed<sup>†</sup>

†SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, we will issue a full refund of the purchase price.  
©2020 SS White All Rights Reserved. SS White® is a registered trademarks of SS White Dental.

**SS WHITE** Dental®

Better Patient Outcomes  
Improved Efficiency  
Faster Practice Growth

1145 Towbin Avenue Lakewood, New Jersey 08701

[www.sswhitedental.com](http://www.sswhitedental.com)