DESCRIPTION

SS White® EndoGuide™ Burs are a unique series of eight burs for non-surgical root canal treatment. EndoGuide™ Burs are specifically designed to maximize efficiency during endodontic exploration and access while conserving healthy tooth structure. The patented¹, conical shaped micro-diameter tip of the EndoGuide™ Burs acts as a self-centering guide to permit straight-line access to canals. EndoGuide™ Burs blade design functions to effectively polish the dentin surface as it cuts, making visual dentin mapping easier and canal identification more precise. Two EndoGuide™ kits address the specific endodontic needs of Anterior/Bicuspid teeth and Molar teeth. Individual EndoGuide™ Burs are also available. (See **Table A**).

EndoGuide™ Anterior/Bicuspid Kit
(For Endodontic Access and Exploration) #18052

Contains all instrumentation to create endodontic access through metal, porcelain and zirconia, featuring SS White® Great White® #2 metal cutting bur, Great White® Z Diamonds along with four EndoGuide™ Burs ideally suited for locating and accessing single root canals in anterior/bicuspid teeth.

EndoGuide™ Molar Kit
(For Endodontic Exploration) #18051

Contains seven EndoGuide™ Burs designed to increase visibility and control during endodontic exploration in molars when locating hidden canals, navigating deeply calcified canals and troughing between canals.

NOTE: Initial access and de-roofing in molars should be performed with either the Great White® #2 bur or Great White® Z Diamonds supplied in the EndoGuide™ Anterior/Bicuspid Kit for Endodontic Access and Exploration, order #18052.

OPERATIONAL INSTRUCTIONS:

- EndoGuide™ Burs are task specific; each bur is designed for a particular area of use during endodontic access and canal exploration.
- While EndoGuide™ Burs are diminutive in size, they are very efficient cutters in both the vertical and lateral directions. To maximize cutting efficiency and avoid bur breakage, it is recommended that the clinician apply only enough pressure to guide the bur while allowing the instrument to progress through the substrate. When used for troughing procedures, a light sweeping motion is suggested.
- Excessive force to the instrument, especially in a lateral direction, can lead to bur damage.

INITIAL ACCESS THROUGH RESTORATIONS/TOOTH SURFACES:

- Identify the material makeup of any existing restorations that must be removed or penetrated to gain access.
- Once identified, select the appropriate instrument for the task. (See **Table A**)

NOTE: SS White® Complete Restoration Removal Kit (Order #18180) contains instrumentation specifically designed for efficient, non-traumatic removal of failed amalgams, composites, PFM crowns, and all-ceramic crowns, including zirconia and lithium disilicate ceramics.

ENDODONTIC ACCESS AND EXPLORATION:

EndoGuide™ Burs are ideal for use when performing micro-endodontics using magnification-enhanced vision via a dental operating microscope or loops. For best product performance, always match the EndoGuide™ Bur to the task to be performed (See **Table A**).

Slow-speed EndoGuide™ Burs are those with SLRA (Surgical Length Right Angle) shanks and XLRA (Extra-Long Right Angle) shanks for use in standard slow-speed right-angle handpieces at slow-speed rpm ranges appropriate for endodontic procedures.

High-speed EndoGuide™ Burs are those with SLRA (Surgical Length Friction Grip) shanks for use in standard high-speed friction grip handpieces at high-speed rpm anges appropriate for endodontic procedures.

PRECAUTIONS:

- Before each use, inspect the bur for damaged blades or signs of excessive wear. Discard worn burs.
- Maintain handpieces in accordance with the manufacturer's instructions.
- Fully seat the bur into the chuck mechanism to eliminate non-concentric operation.
- Do not exceed manufacturer's suggested speed recommendations.
- Clean and sterilize burs in accordance with the instructions provided by the manufacturer of your sterilization unit. Do not immerse carbide burs in cold sterilizing solutions or other strong oxidizing agents.

WARNINGS:

SS White® EndoGuide™ Burs may not be appropriate for all patients under all circumstances, so professional judgment and discretion should be exercised at all times. SS White® EndoGuide™ Burs should be used only by licensed dental professionals with appropriate levels of clinical experience and training. Caution and care should be used at all times. Treatment options and techniques outlined in this insert are suggestions only and may not be appropriate for all patients. Dentists must apply appropriate levels of force when using SS White® EndoGuide™ Burs. Excessive force can cause premature wear to the instruments. Use your professional judgment.

LIMITATION OF LIABILITY

SS White® disclaims all liability and damages of any kind resulting from intentional or negligent conduct or otherwise that arise out of the use or misuse of any SS White® instruments or these instructions.

1145 Towbin Avenue, Lakewood, NJ 08701 • 732-905-1100 / 800-535-2877
www.sswhiteburs.com
1. U.S. Patent No. 6,257,889 B1
SS White® and Great White® are registered trademarks of SS White Burs, Inc. All rights reserved.
SS White® is a registered ISO 13485 Company. © 2010 SS White Burs, Inc.
Ver. 2 - 7/2011

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Ver. 1.1 - 7/2011

Use the following reference guide (Table A.) when matching EndoGuide™ Burs to the task to be performed

Kit and Bur Reference Guide: SS White® EndoGuide™ Precision Micro Endodontic Burs

Bur#	Bur Photograph (Actual Size)		Head Length (mm)	Reorder Pkg. Size	SS White Order #	Included in: Anterior/Bicuspid Kit#18052	Included in: Molar Kit Kit#18051	Recommended Applications
EG 1 A		SLFG* Length = 27mm	2.5	5-pk	15066	YES	NO	• Initial access in non-restored anterior and bicuspid teeth
EG 1		SLRA* Length = 27mm	3.5	5-pk	14780	YES	YES	• Deep troughing • Deeper orifice enlargement and calcified canals
EG 2		SLRA* Length = 27mm	2.5	5-pk	14779	YES	YES	• De-roofing pulp chamber • Deep orifice enlargement and calcified canals
EG 3		SLFG* Length = 27mm	1.5	5-pk	15023	YES	YES	• Initial access for small incisors • Troughing and navigating calcified canals
EG 4		SLFG* Length = 29mm	3.5	5-pk	15029	NO	YES	• Deep troughing • Navigating super ovoid and calcified canals
EG 5		XLRA* Length = 34mm	1.5	5-pk	14794	NO	YES	• Deep troughing
EG 6		XLRA* Length = 34mm	2.5	5-pk	14795	NO	YES	• Retrieving separated instruments
EG 7		XLRA* Length = 34mm	1.5	5-pk	14786	NO	YES	• Navigating deeply calcified canals
SS White Great White® Carbides and Diamonds:								
GW2		FG Length = 19mm	4.0	10-pk 100-pk	15062 13095	YES	NO	• Metal crown access or removal • Amalgam restoration and decay removal • Molar initial access and de-roofing
GWZ		FG Length = 22mm	8.0	5-pk	18161	YES	NO	• Ceramic/zirconia crown removal • Molar initial access and de-roofing
GWZ		FG Length = 19mm	1.8	5-pk	18163	YES	NO	• Initial access through ceramic/zirconia crowns • Molar initial access and de-roofing

*NOTE: SLRA (Surgical Length Right Angle) and XLRA (Extra-Long Right Angle) shanks are for use in right angle slow-speed handpieces; SLFG (Surgical Length Friction Grip) shanks are for use in friction grip high-speed handpieces.

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