



ORAL SURGERY BURS























*The Fastest Cutting Burs
Engineered for maximum strength and durability
to save you time and money*



simply better dentistry™

WHITE®

ORAL SURGERY BUR SHAPE SELECTION (SOLD IN 10 PACKS)

			Shank #1		Shank #2	Shank #3	Shank #4	Shank #5	Shank #6			
Round	SIZE	HEAD DIA.	SSW#	Sterile Package*	SSW#	SSW#	SSW#	STRYKER#	SSW#	STRYKER#	SSW#	STRYKER#
	1	1.0	14823									
	2	1.2	14826									
	4	1.4	14832									
	6	1.8	14838	34838	30020	30035	30050	2296-100-418	30060	2296-101-418	30000	277-10-418
	8	2.3	14844	34844	30021	30036	30051	2296-100-423	30061	2296-101-423	30001	277-10-423
	10	2.7	14846	34846								
Cross Cut Fissure												
	557	1.0	14871									
	558	1.2	14874									
	560	1.6	14880				30052				30002	
	1557	1.0	14898									
	1558	1.2	14901									
	700	1.2	14883									
	701	1.2	14886	34916	30022		30053	2296-100-212	30062			
	701L	1.2	14910									
	702	1.6	14889	34889	30023	30037	30054	2296-100-216	30063	2296-101-216	30003	277-10-216
	702L	1.6	14890		30024		30055	2296-100-316	30064	2296-101-316	30004	
	1702	1.6	14916		30029		30048		30069			
	1702L	1.6	14917		30030		30049		30070			
	703	2.1	14892	34892	30025	30038	30056	2296-100-221	30065	2296-101-221	30005	277-10-221
	703L	2.1	14893		30026	30039	30057	2296-100-321	30066	2296-101-321	30006	
	1703	2.1	14914	34914	30027	30040	30058		30067		30007	
	1703L	2.1	14915	34915	30028		30059	2296-100-421	30068	2296-101-421	30008	

(SOLD IN 5 PACKS)

Surgical Bur (SL) Shank Length - Round				Surgical Bur (SL) Shank Length - Cross Cut				Surgical Bur (SL) Shank Length - Cross Cut			
SIZE	HEAD DIA.	SLRA	SLFG	SIZE	HEAD DIA.	SLRA	SLFG	SIZE	HEAD DIA.	SLRA	SLFG
1	1.0		14001-5	557	1.0	14157-5	14057-5	701	1.2		14071-5
2	1.2	14112-5	14002-5	558	1.2	14158-5	14058-5	702	1.6	14172-5	14072-5
4	1.4	14114-5	14004-5	1557	1.0	14159-5	14059-5	703	2.1		14073-5
6	1.8	14116-5	14006-5	1558	1.2		14060-5	1702	1.6	14174-5	14074-5
8	2.3	14118-5	14008-5	700	1.0		14070-5	1703	2.1	14175-5	14075-5


SS White® Oral Surgery Burs Are...

- Manufactured to exceed ISO specifications for quality and consistency
- Packaged sterile for safety and convenience
- Compatible with all popular hand pieces including Stryker®, Hall® and MicroAire®

A recent study confirms that SS White® burs are 84.67% stronger than competitor burs.* SS White® goes above and beyond ISO and ADA testing requirements. The neck can withstand 16.41 (lbs) of force, which is 40% stronger than competitors.*

*Internal testing data available upon request

Third Molar Extractions



Reduced Clogging

SS White®'s unique blade geometry allows for more effective bone and tooth removal discharge creating a more efficient instrument

More Efficient Cutting

Faster cutting by design, SS White® Oral Surgery Burs save time for all surgical procedures





Reduced Bur Breakage

Internal testing shows that SS White® burs consistently outperform all competitive burs in neck strength, significantly reducing bur breakage

Packaged Sterile

For your convenience and safety, SS White® Oral Surgery burs shank #1-6* come in a single sterile package

For Precious and Non-Precious Crown and Bridge

GREAT WHITE®				
SIZE/FG	GW2(1558)	GW2S	GWH342	GWH342L
10 Pack	15062	15064	14919	14920
Bulk Pack	13095	13096		

Innovative Design...

SS White® continually looks to improve our Oral Surgery Burs. An example of this is the 1703, which combines the sculpting ability of a round bur with the high performance of a 703 sectioning bur. Now, you can have the best of both worlds in one bur.

At SS White® we are constantly striving to exceed your expectations.



Oral Surgery Bur Shank Length #1-6

Shank #1 - 44.5mm

Fits Hall® & Most Standard Hand Pieces for Medium Bur Guard



Shank #2 - 51mm

Fits Hall® & MicroAire® Hand Pieces with Medium Bur Guard



Shank #3 - 65mm

Fits Hall® & MicroAire® Hand Pieces with Long Bur Guard



Shank #4 - 44.5mm

Fits Stryker® 100K & Impaction Handpiece



Shank #5 - 59mm

Fits Stryker® 100K & Impaction Handpiece



Shank #6 - 44.5mm

Fits Stryker® 50K & Osteotome Drill J-Notch



Surgical Bur (SL) Shank Length

Friction Grip Surgical Length



.98"
25mm LOA


Right Angle Surgical Length



1.023"
26MM LOA

simply better dentistry™
SS WHITE®

Oral Surgery Burs Brochure - Rev. 2



6 4 6 0 7

800.535.2877 | www.sswhitedental.com
1145 Towbin Avenue Lakewood, New Jersey 08701