



ORTHODONTIC PRODUCTS

SS White Orthodontic Carbides and Diamonds



simply better dentistry™


SS WHITE®

Orthodontic Finishing with Consistent and Efficient Results

SS White Orthodontic Carbides and Diamonds enable faster, more predictable orthodontic debonding and finishing procedures, creating better patient outcomes with beautiful smiles.

Reliant™ Orthodontic Finishing Kit

High-quality carbide debonding burs allow predictable removal of cement without damage to the enamel, with concave shapes convenient for convex buccal surfaces.



Kit Includes (1 of each):
 118S-18 Blade Carbide Finishing Bur #15122-5
 118L-18 Blade Carbide Finishing Bur #15123-5
 218-18 Blade Carbide Finishing Bur #15124-5
 129- Diamond (868-018EF) #91129-5
 815SL- 18 Blade Finishing Bur #15125-5
 Jazz Polisher Universal Flame #89031
 Kit Reorder No. #15084

JAZZ
POLISHERS

Jazz Supreme


Polish enamel after debonding. Made with specially formulated rubber matrix providing an excellent polish with minimal pressure.



Jazz® Polishers Supreme Composite 1-Step Universal Polishing System 4 Pc. Kit in Autoclavable Bur Block
 Kit Includes (1 of each):
 #89031 Universal Flame
 #89032 Universal Knife Edge
 #89033 Universal Cup
 #89034 Universal Small Flame
 Kit Reorder No. #89035

TRIMMING & FINISHING BURS

Perfect for removing bonding adhesive and cement after bracket removal, creating ultra-smooth, pre-polishes.



Perfect for removing bonding adhesive and cement after bracket removal, creating ultra-smooth, pre-polishes.

FOOTBALL								FOA1		BULLET					
	12 BLADE	15106-5	15095-5	15408-5			12 BLADE		16004-5		12 BLADE	15801-5	15802-5	15092-5	14333-5
	20 BLADE				15404-5	15406-5					30 BLADE				
	30 BLADE					15409-5					BULK 100			13626	
BULK 100	13623	13624	13625				14316-5								

FLAME						INVERTED TAPER			SAFE END			
	12 BLADE	15091-5	15094-5	15108-5	15274-5		12 BLADE	15303-5		15304-5	30 BLADE	14361-5
	20 BLADE											
					15275-5							

ROUND END TAPER					TAPER						
	12 BLADE	15675-5	15677-5			12 BLADE	15705-5	15713-5	15642-5	15664-5	
	20 BLADE			15676-5		15678-5					14082-5
											14342-5

FLAT END TAPER							CFT				
	12 BLADE	15204-5	1520-5	15206-5	15275-5	15381-5		FG	16001-5	16002-5	16003-5
	20 BLADE										

SURGICAL LENGTH FINISHING					SAFE END SERIES							
	12 BLADE	16902-5	16900-5	14087-5		16901-5	10 BLADE	16006-5	16007-5	16008-5	16909-5	16010-5
							20 BLADE	16030-5	16031-5	16032-5	16033-5	16034-5

Delivers minimal vibration to maximize performance and procedural comfort for you and your patient.

											
	4	6	35	330	56	57	557	1157	1557	169L	699L
FG	14970	14973	14985	15001	14991	14994	15006	17708	15021	16884	16890
FG (BULK)	13004	13006	13035	13330	13056	13057	13557	13157	13559	13169	13699

GREAT WHITE[®] GOLD SERIES

METAL CUTTING - CAVITY PREPARATION CARBIDE BURS

	
	GW2
FG	15062
FG (BULK)	13095

GREAT WHITE[®] ULTRA

METAL APPLIANCE CUTTING AND TRIMMING BURS

		
	GWU856-020	GWU856-023
FG	18189-5	18190-5

PIRANHA[®] SINGLE-PATIENT-USE DIAMOND INSTRUMENTS

Easily smooth and contour contact areas after stripping, with wide selection of short shanks and shorter-than-short shanks for improved access.

							
SHAPE	368-016	392-016	858-008	858-014	859-010	859-014	862-012
DIAMETER (mm)	1.6	1.6	0.8	1.4	1.0	1.4	1.2
HEAD LENGTH (mm)	3.4	4.5	8.0	8.0	10.0	10.0	8.0
TIP DIAMETER (mm)	-	-	0.3	0.5	0.6	0.25	0.35
● Super Coarse				858-014SC			862-012SC
● Coarse	368-016C			858-014C	859-010C	859-014C	862-012C
● Medium	368-016M	392-016M	858-008M	858-014M	859-010M	859-014M	862-012M
● Fine	368-016F	392-016F	858-008F	858-014F	859-010F	859-014F	862-012F
● Very Fine	368-016VF	392-016VF	858-008VF	858-014VF		859-014VF	862-012VF
S (Short Shank)				858-014C-S		859-014C-S	862-012C-S
SS (Shorter than Short Shank)							862-012C-SS
							862-012MSS

Shank Lengths

from 19 to 24 mm	from 16 to 21 mm	16 mm
STANDARD FG	S = SHORT SHANK	SS = SHORTER-THAN-SHORT SHANK

Grit Size Designation

● Super Coarse (SC)	● Fine (F)
● Coarse (C)	● Very Fine (VF)
● Medium (M)	○ Ultra Fine (UF)

Micron (μm)

● 149÷177 (SC)	● 38÷45 (F)
● 105÷149 (C)	● 30÷38 (VF)
● 62÷125 (M)	○ 8÷12 (UF)

REVELATION[®] DIAMOND INSTRUMENTS

Premium Quality - Multi-Use Diamonds




3-layer dense, uniform diamond matrix keeps greater number of diamond particles in direct contact with tooth surface for smooth feel and consistently fast cutting, providing less vibration and greater comfort for the patient.

						
FG	7901-009	7404-014	Safe End	SE3	SE6	SE9
Shape	Finishing					
HEAD DIAMETER (MM)	0.9	1.4		1.0	1.0	1.6
HEAD LENGTH (mm)	3.7	2.6		3.0	6.0	9.0
TIP DIAMETER (mm)	-	-		0.9	0.9	0.6
● Fine	91154-5	91156-5		91148-5	91150-5	91152-5
● Very Fine	91155-5	91157-5		91149-5	91151-5	91153-5




EXPRESS LINE DIAMOND DISCS

Uniform spacing between the diamond aggregates helps ensure immediate removal of debris and less clogging for orthodontic stripping.


(F) FLEX Mounted Diamond Discs



(S) SUPERFLEX Mounted Diamond Discs




Standard Straight Handpiece



GREAT WHITE[®] LAB SERIES CARBIDE BURS

Abrasion resistance reduces surface heat and vibration for a cooler, more consistent surface finish, and Zirconium Nitride coating provides increased performance and durability.



Ortho Product Brochure



6 4 6 6 3

800.535.2877 | www.sswhitedental.com
1145 Towbin Avenue Lakewood, New Jersey 08701