
World Wide DIAMOND Instruments



SS WHITE®

The Brand Trusted by Clinicians Worldwide
Our Instruments, Your Artistry



World Wide DIAMOND™ Instruments

Grit Size Designation

- Coarse (C) ● Very Fine (VF)
- Medium (M) ○ Ultra Fine (UF)
- Fine (F)

Micron (µm)

- 105-149 (C) ● 30-38 (VF)
- 62-125 (M) ○ 8-12 (UF)
- 38-45 (F)

Shank Length

from 19 to 24 mm
STANDARD FG

All available
in 10 pack
blister packs

Round

SHAPE	W801-008	W801-009	W801-010	W801-012	W801-014	W801L-014	W801-016	W801L-016	W801-018	W801L-018	W801-021	W801L-021	W801-023	W801-035	W802-009
DIAMETER (mm)	0.8	0.9	1.0	1.2	1.4	1.4	1.6	1.6	1.8	1.8	2.1	2.1	2.3	3.5	0.8
HEAD LENGTH (mm)	0.7	0.8	0.9	1.1	1.3	1.3	1.5	1.5	1.7	1.7	2.0	1.9	2.2	3.4	2.2
TIP DIAMETER (mm)	-	-	-	-	-	-	-	-	-	-	-	-	-	-	0.9
● Coarse				W801-012C	W801-014C			W801L-016C	W801-018C	W801L-018C	W801-021C				
● Medium	W801-008M	W801-009M	W801-010M	W801-012M	W801-014M	W801L-014M	W801-016M	W801L-016M	W801-018M	W801L-018M	W801-021M	W801L-021M	W801-023M	W801-035M	W802-009M
● Fine									W801-018F						
● Very Fine					W801-014VF								W801-023VF		

Round

Inverted Cone

SHAPE	W802-010	W802-012	W802-014	W802-016	W802-018	W805-010	W805-012	W805-014	W805-016	W805-018	W805-021	W807-014	W807-016	W807-018
DIAMETER (mm)	1.0	1.0	1.0	1.2	1.2	0.8	0.9	1.1	1.2	1.3	1.5	1.2	1.4	1.5
HEAD LENGTH (mm)	2.5	2.6	3.0	3.0	3.6	1.1	1.2	1.4	1.6	1.8	2.1	3.0	4.0	5.0
TIP DIAMETER (mm)	1.0	1.2	1.4	1.6	1.8	1.0	1.2	1.4	1.6	1.8	2.1	1.4	1.6	1.8
● Coarse			W802-014C					W805-014C		W805-018C		W807-014C		
● Medium	W802-010M	W802-012M	W802-014M	W802-016M	W802-018M	W805-010M	W805-012M	W805-014M	W805-016M	W805-018M	W805-021M	W807-014M	W807-016M	W807-018M

Double Inverted Cone

Wheel

Pear

SHAPE	W806-012	W806-014	W806-016	W815-016	W818-045	W819-035	W909-035	W909-042	W830-008	W830-010	W830-012	W830L-014	W830-016	W830-021	W830L-021
DIAMETER (mm)	0.8	1.0	0.9	1.6	4.5	3.5	3.5	4.2	0.7	0.8	0.9	1.1	1.3	1.8	1.2
HEAD LENGTH (mm)	2.9	3.0	3.0	0.5	0.5	0.9	1.5	1.2	2.5	3.0	3.0	5.0	3.0	5.0	5.0
TIP DIAMETER (mm)	1.2	1.4	1.6	-	-	-	-	-	0.8	1.0	1.2	1.4	1.6	2.1	2.1
● Coarse										W830-010C					
● Medium	W806-012M	W806-014M	W806-016M	W815-016M	W818-045M	W819-035M	W909-035M	W909-042M	W830-008M	W830-010M	W830-012	W830L-014M	W830-016M	W830-021M	W830L-021
● Fine												W830L-014F			

Barrel

Tri-Wheel Depth Cutter

Flat-End Cylinder

SHAPE	W811-031	W811-047	W834-016	W834-021	W835-008	W835-010	W835-012	W835-014	W835-016	W836-012	W836-014	W837-012	W837-014	W837-016	W837-018
DIAMETER (mm)	3.1	4.7	1.6	2.1	0.8	1.0	1.2	1.4	1.6	1.2	1.4	1.2	1.4	1.6	1.8
HEAD LENGTH (mm)	4.0	7.0	6.0	6.0	3.0	4.0	4.0	4.0	4.0	6.0	6.0	8.0	8.0	8.0	8.0
TIP DIAMETER (mm)	-	-	-	-	0.8	1.0	1.2	1.4	1.6	1.2	1.4	1.2	1.4	1.6	1.8
● Coarse							W835-012C			W836-012C	W836-014C	W837-012C	W837-014C	W837-016C	W837-018C
● Medium	W811-031M	W811-047M	W834-016M	W834-021	W835-008M	W835-010M	W835-012M	W835-014M	W835-016M	W836-012M	W836-014M	W837-012M	W837-014M	W837-016M	W837-018M

Round-End Cylinder

Christmas Tree

Flat-End Taper

SHAPE	W838-010	W838-012	W880-012	W880-014	W880-021	W881-010	W881-012	W881-014	W881-016	W852-010	W852-012	W844-012	W844-014	W845-014	W846-014
DIAMETER (mm)	1.0	1.2	1.2	1.4	2.1	1.0	1.2	1.4	1.6	1.0	1.2	1.2	1.4	1.4	1.4
HEAD LENGTH (mm)	4.0	4.0	6.0	6.0	6.0	8.0	8.0	8.0	8.0	6.0	6.0	6.0	6.0	4.0	6.0
TIP DIAMETER (mm)	1.0	1.2	1.2	1.4	2.1	1.0	1.2	1.4	1.6	0.5	0.5	0.9	0.9	0.9	0.7
● Coarse		W838-012C	W880-012C					W881-014C		W852-010C					
● Medium	W838-010M	W838-012M	W880-012M	W880-014M	W880-021M	W881-010M	W881-012M	W881-014M	W881-016M	W852-010M	W852-012M	W844-012M	W844-014M	W845-014M	W846-014M
● Very Fine												W852-012VF			



SS WHITE®

732-905-1100 | www.sswwhiteburs.com

1145 Towbin Avenue • Lakewood, NJ 08701

World Wide Diamond Instruments Catalog



6 4 6 6 2