

SS White - Guía de aplicación para las fresas CAD CAM para molinos de Roland

Para garantizar un uso eficiente de las fresas SS White CAD CAM, consulte nuestra guía de aplicación para seleccionar las herramientas correctas para su uso en su molino dental de Roland según el modelo y el material a fresar. Consulte las instrucciones del fabricante Roland para asegurar el uso adecuado de las fresas.

Igualmente importante es seleccionar el conjunto de herramientas correcto dentro de la clave de producto (software). Las herramientas SS White están diseñadas para garantizar la funcionalidad plug and play cuando se utiliza con las especificaciones predeterminadas de fábrica (osición, velocidad y dimensiones). Si hay alguna incertidumbre sobre la coincidencia dimensional entre las herramientas, mida y compare las longitudes de las puntas antes de su uso. Si no lo hace, podría dañar la herramienta y la máquina.

Disponible a través de: **SSWHITE Dental®**

SSW Numero de Artículo	Tipo	Diámetro de Cabeza	Diámetro del Vástago	Diseño de Cortador	Máquina	Material
 SSWR0L-034050-20-02	DuraDi	.3mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
	DuraDi	.6mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
	DuraDi	.8mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
	DuraDi	1mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
	DuraDi	2mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
 SSWR0L-034050-20-01	Lazer Sharp	.3mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia
	Lazer Sharp	.6mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia
	Lazer Sharp	.8mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia
	Lazer Sharp	1mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia
	Lazer Sharp	2mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia