

SS White - Roland Cad Cam Milling Bur Application Guide

To ensure seamless use of the SS White Cad Cam milling burs, please consult our application guide to select the correct tools for use in your Roland dental mill according to the model and material to be milled. Consult the Roland dental mill manufacturer instructions to insure proper use of the milling burs.

Equally important is selecting the correct tool set within your software package. SS White tools are designed to ensure plug and play functionality when used with factory default specifications (feed, speed, and dimensions). If there is any uncertainty about the dimensional match between the tool and tool library, measure and compare tip lengths prior to use. Failure to do so could damage the tool, substrate and or machine.

Available Through:

SS WHITE Dental®

SSW Item #	Type	Head Diam.	Shank Diam.	Cutter Design	Machine	Material
						
SSWR0L-034050-20-02	DuraDi	.3mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
SSWR0L-064050-20-02	DuraDi	.6mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
SSWR0L-084050-20-02	DuraDi	.8mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
SSWR0L-104050-20-02	DuraDi	1mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
SSWR0L-204050-20-02	DuraDi	2mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia, Some Composites
						
SSWR0L-034050-20-01	Lazer Sharp	.3mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia
SSWR0L-064050-20-01	Lazer Sharp	.6mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia
SSWR0L-084050-20-01	Lazer Sharp	.8mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia
SSWR0L-104050-20-01	Lazer Sharp	1mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia
SSWR0L-204050-20-01	Lazer Sharp	2mm	4mm	2 Flute Ball End	Roland DWX-4,DWX-50, 51D, 52DC	PMMA,Wax,Zirconia