

INDIRECT RESTORATIVE MATERIALS

ZIRCONIA - Monolithic or Layered

Bruxzir

Crown Removal

Removal	Grit	
Great White®Z 856-018M	M	
Great White®Z 850-018M	M	

Endo Access

Great White®Z 801-016M	M	
Great White®Z 801-018M	M	

Preparation	Grit	
Diamond 856-018M (Axial Wall Reduction)	M	
Diamond 858-014M (Interproximal Reduction)	M	
Diamond 379-023M (Occlusal/Lingual Reduction)	M	
Diamond 909-035M (Incisal Reduction)	M	
Carbide Alternative		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-023 (Axial Wall Reduction)		
Great White® Ultra Carbide 379-023 (Occlusal/Lingual Reduction)		

Adjust/Pre-Polish	Grit	
Adjust		
Diamond 868-016F	F	
Diamond 868-023F (Occlusal/Lingual)		
Diamond 850-016F (Margin Refinement)	F	
Diamond 859-016F (Interproximal)	F	
Pre-Polish:		
Diamond 868-018VF	VF	
Diamond 868-023VF (Occlusal/Lingual)		
Diamond 850-015VF (Margin Refinement)	VF	
Diamond 859-016VF (Margin Refinement)	VF	

Polish	Grit	
Jazz® polisher Medium Cup ZA2S	M	
Jazz® polisher Fine Cup ZA2S	F	
Jazz® polisher Medium Flame ZA2S	M	
Jazz® polisher Fine Flame ZA2S	F	

PORCELAIN/MILLED CERAMICS

Feldspathic leucite reinforced porcelain (Empress) - Lithium disilicate porcelain crowns (Emax) - LAVA or Procera

Crown Removal

Removal	Grit	
Great White®Z 856-018M	M	
Great White® Z 850-018M	M	

Endo Access

Great White®Z 801-016M	M	
Great White®Z 801-018M	M	

Preparation	Grit	
Diamond 878K-016M	M	
Diamond 878K-018M	M	
Diamond 878K-021M (Axial Wall Reduction)	M	
Diamond 878K-014M (Axial Wall Reduction)	M	
Diamond 368-023M (Occlusal/Lingual Reduction)	M	
Diamond 909-035M (Occlusal/Lingual Reduction)	M	

Adjust/Pre-Polish	Grit	
Diamond 368-016F	F	
Diamond 368-023F (Occlusal/Lingual)	F	
Diamond 878K-016F	F	
Diamond 878K-018F		
Diamond 878K-021F (Marginal Refinement)		
Diamond 878K-014F (Interproximal)	F	
Pre-Polish:		
Diamond 368-016VF	VF	
Diamond 368-023VF (Occlusal/Lingual)		

Polish	Grit	
Jazz® polisher Medium Cup P3S	M	
Jazz® polisher Fine Cup P3S	F	
Jazz® polisher Medium Flame P3S	M	
Jazz® polisher Fine Flame P3S	F	

PORCELAIN

(Fused-to Metal)

Ceramic Portion

Removal	Grit	
Porcelain		
Diamond 856-016M	M	
Diamond 856-018M	M	
Metal		
Great White® #2		
Great White® #1557		
Decay Removal		
FG #1557		
FG #1558		
FG #330		
FG #245		

Metal Portion

Preparation	Grit	
Diamond 856-018M	M	
Diamond 850-018M (Axial Wall Reduction)		
Diamond 858-014M (Interproximal Reduction)	M	
379-023M (Occlusal/Lingual Reduction)	M	
909-035M (Incisal Reduction)	M	
Carbide Alternative		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-023 (Axial Wall Reduction)		
Great White® Ultra Carbide 379-023 (occlusal surface)		

Adjust/Pre-Polish	Grit	
Adjust		
Diamond 379-018F	F	
Diamond 379-023F (Occlusal/ Lingual)		
856-018F	F	
850-018F (Margin Refinement)		
859-016F (Interproximal)	F	
Pre-Polish:		
379-018VF	VF	
379-023VF (Occlusal/ Lingual)		
856-015VF (Margin Refinement)	VF	
850-015VF (Margin Refinement)	VF	
859-016VF (Interproximal)	VF	

Polish	Grit	
Jazz® polisher Medium Cup P3S	M	
Jazz® polisher Fine Cup P3S	F	
Jazz® polisher Medium Flame P3S	M	
Jazz® polisher Fine Flame P3S	F	
Jazz® polisher Coarse Cup P3S	C	
Jazz® polisher Coarse Cup P3S	C	

GOLD ALLOYS
(Nobel)

Removal		
Porcelain		
Diamond 856-016M	M	
Diamond 856-018M	M	
Metal		
Great White® #2		
Great White®#1557		
Decay Removal		
FG #1557		
FG #1558		
FG #330		
FG #245		

Preparation		
Diamond 856-018M (Axial Wall Reduction)	M	
858-014M (Interproximal Reduction)	M	
379-023M (Occlusal/ Lingual Reduction)	M	
909-035M (Incisal Reduction)	M	
Carbide Alternative		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-023 (Axial Wall Reduction)		
GWU 379-023 (Occlusal/ Lingual Reduction)		

Adjust/Pre-Polish		
Adjust		
Diamond 868-016F	F	
Diamond 868-023F (Occlusal/ Lingual)		
Diamond 850-016F (Margin Refinement)	F	
Diamond 859-016F (Interproximal)	F	
Pre-Polish:		
Diamond 868-018VF	VF	
Diamond 868-023VF (Occlusal/ Lingual)		
Diamond 850-015VF (Margin Refinement)	VF	
859-016VF (Interproximal)	VF	

Polish		
Jazz® polisher Medium Cup P3S	M	
Jazz® polisher Fine Cup P3S	F	
Jazz® polisher Medium Flame P3S	M	
Jazz® polisher Fine Flame P3S	F	

NICKEL OR
COBALT-CHROME
(Base-Metal)
Alloys
or
Titanium
Abutments

Crown Removal		
Porcelain		
Diamond 856-016M	M	
Diamond 856-018M	M	
Metal		
Great White® #2		
Great White® #1557		
Decay Removal		
FG #1557		
FG #1558		
FG #330		
FG #245		

Preparation		
Diamond 856-018M (Axial Wall Reduction)	M	
858-014M (Interproximal Reduction)	M	
379-023M (Occlusal/ Lingual Reduction)	M	
909-035M (Incisal Reduction)	M	
Carbide Alternative		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-023 (Axial Wall Reduction)		
GWU 379-023 (Occlusal/ Lingual Reduction)		

Adjust/Pre-Polish		
Adjust		
Diamond 868-016F	F	
Diamond 868-023F (Occlusal/ Lingual)		
Diamond 850-016F (Margin Refinement)	F	
Diamond 859-016F (Interproximal)	F	
Pre-Polish:		
Diamond 868-018VF	VF	
Diamond 868-023VF (Occlusal/ Lingual)		
Diamond 850-015VF (Margin Refinement)	VF	
859-016VF (Interproximal)	VF	

Polish		
Jazz® polisher Medium Cup P3S	M	
Jazz® polisher Fine Cup P3S	F	
Jazz® polisher Medium Flame P3S	M	
Jazz® polisher Fine Flame P3S	F	

DIRECT RESTORATIVE MATERIALS

AMALGAM

Restorative Removal	Grit	Preparation	Grit	Adjust/Contour	Grit	Polish	Grit
Amalgam Removal: Great White® #1	M	Diamond 835-010M	M	Adjust: Diamond 868-016F Diamond 868-023F (Occlusal)	F	Jazz® polisher Medium Cup P3S	M
Great White® #2		Diamond 835-012M	M	Pre-Polish Diamond 868-018VF Diamond 868-023VF (Occlusal)	VF	Jazz® polisher Fine Cup P3S	F
Diamond 813-014M		Diamond 845-010M	M	Carbide Alternative: Adjust: 7404 7406 7408	VF	Jazz® polisher Medium Flame P3S	M
Caries Removal: RA #4		Diamond 845-012M	M	Pre-Polish: 8404 8406 8408		Jazz® polisher Fine Flame P3S	F
RA #6		Carbide Alternative #556				Jazz® polisher Coarse Cup P3S	C
FG #4		#557				Jazz® polisher Coarse Cup P3S	C
FG #6		#558					
FG #330							
FG #1557							
Diamond 830-012M	M						

DIRECT AND INDIRECT RESTORATIVE MATERIALS

COMPOSITE
GLASS
IONOMER
restoration

RESIN
IONOMER
restoration

Restorative Removal	Grit	Preparation	Grit	Adjust/Contour	Grit	Polish	Grit
Amalgam Removal: Great White® #1	M	Diamond 830-010M	M	Adjust: Diamond 868-016F Diamond 868-023F (Occlusal)	F	Jazz® polisher Fine Cup C1S	F
Great White® #2		Diamond 830-012M	M	Pre-Polish Diamond 868-018VF Diamond 868-023VF (Occlusal)	VF	Jazz® polisher Medium Flame C1S	M
Diamond 813-014M		Diamond RA 801-010M	M	Carbide Alternative: Adjust: 7404 7406 7408	VF		
Caries Removal: RA #4		Diamond RA 801-012M	M	Pre-Polish: 8404 8406 8408			
RA #6		Carbide Alternative #330					
FG #4		#245					
FG #6		#1557					
FG #330		#558					
FG #1557		Fissurotomy® Original					
Diamond 830-012M	M						

#64669

Better Patient Outcomes
Improved Efficiency
Faster Practice Growth

SS WHITE Dental®

1145 Towbin Avenue • Lakewood, NJ
sswhitedental.com

Dentists Materials
Reference



Better Patient Outcomes
Improved Efficiency
Faster Practice Growth

